

Non-invasive Interactive Neurostimulation (InterX®) elicits significantly greater physiological response than TENS: *Lymphocyte Metabolism and Cytokine production*

Introduction:

Non-invasive interactive neurostimulation (InterX) has shown to be clinically effective for a range of acute and chronic conditions^{1,2,3}and in experimental pain research⁴. InterX also demonstrated reduced inflammation in patients following ankle surgery². While mechanisms for pain relief for this type of therapy are well understood^{5,6}, it has previously been demonstrated that transcutaneous electrical stimulation does not affect the inflammatory response⁷.

In a previous study, the effects of InterX pre- and post-treatment were examined. It was demonstrated that lymphocytes post treatment may be better able to respond to stimuli when necessary, which could result in enhanced healing ability post injury. The genomic profile seen after just 20 minutes is indicative of an acute response that would infer improved ability of these cells to respond to injury types of stimuli. We extended this study to include volunteers treated with a TENS device to compare physiological response following a different form of electrical stimulation.


Hypothesis:

In this study it is hypothesized that significantly greater physiological response will occur following InterX treatment than following TENS. This is suggested in view of the higher amplitude (approximately 3 - 5 times higher) and current density (approximately 90 times higher) output of the InterX device.

Method:

Treatment protocol: A multidisciplinary approach was used to examine the effect of transdermal interactive neurostimulation using the InterX 5002 device compared to the Biomed 2000XL TENS device on lymphocyte metabolic function and cytokine production. Blood was drawn from 8 healthy adults (4M/4F) prior to and 20 minutes following a treatment session. Treatment applications consisted of 10 minutes of treatment on the lateral elbow of the arm from which blood was drawn and 10 minutes over the corresponding spine root. The TENS device had two channels so that the direct area around the elbow, as well as the corresponding spine root could be treated simultaneously to follow the InterX protocol (Fig 1). The manufacturer recommendations to increase the intensity to a level just below that which causes muscle contraction were followed. 100Hz frequency was used with a pulse width of 150 microseconds. These parameters matched the InterX protocol as closely as possible. Using a higher frequency (250Hz) caused greater muscle contraction which thus meant using a lower amplitude so in the interests of comparison, frequency was sacrificed in favor of amplitude. The InterX device has small, closely spaced electrodes which allow for the use of high amplitude stimulation without causing muscle contraction. As TENS electrodes need to be larger and are spaced further apart, muscle contraction occurs at much lower amplitude.

Fig 1. InterX and TENS treatment areas



Blood samples were obtained from consented adult volunteers. Non red blood cells (mitochondrial containing cells) were isolated and placed in the chamber of a high resolution respirometry machine (Oroboros, Innsbruck, Austria). All experiments were performed in duplicate. Respiration of the cells in the presence of Glutamate (a favored substrate of lymphocytes) and 2, 4-Dinitrophenol (mitochondrial uncoupling for maximum respiration) was recorded. In order to account for differing numbers of the cell in the chamber, a BCA protein assay was performed to determine total protein to which the respiration data was normalized.

Results:

The TENS treatment when compared to pre-treatment values, slightly but not significantly increased the maximal respiration in the uncoupled state and there was no discernible difference due to Glutamate over control ($p>.05$) (Fig. 2)

The InterX treatment however elicited significantly ($p<0.05$ by ANOVA) greater respiration both in response to Glutamate and DNP compared with pre-treatment control and compared with TENS.


FIG.3 shows the percentage change in lymphocyte cytokine content, expressed as percent change from pre-treatment to 20 minutes post-treatment. N=4 for InterX, N=4 for TENS. Two lymphocyte extractions were performed on each subject's blood sample and the assays were performed in triplicate.

IL-1 β has been shown to increase mitochondrial oxidative metabolism in pancreatic islet cells in vitro, concomitant to the release of insulin and increased glucose oxidation⁸. IL-6 and TNF- α have both been implicated in liver hypermetabolism, eliciting significantly increased mitochondrial oxidative metabolism⁹. Increased IL-1 β levels, both in arterial wall and systematic circulation, accelerated vascular neointimal formation after injury, and suppressed matrix-metalloproteinase-2 expression (and therefore further vessel wall damage) 6 days following injury¹⁰. In mice lacking the gene for IL-6, bone fracture callus formation, mineralization and remodelling are significantly delayed, although by 6 weeks, healing had progressed to the same level as in wildtype mice. This indicates a role for IL-6 in early fracture healing¹¹. IL-6 has also been shown to play a pivotal role in the healing of CNS following trauma. The effect of IL-6 appear to be most evident in revascularization and prevention of blood-brain barrier breakdown¹².

Joseph F. Clark, ATC, PhD¹ Paul J Magee (Neuro Resource Group, Plano, TX). Gail Pyne-Geithman, D Phil² Department of Neurology,¹ Department of Neurosurgery², University of Cincinnati, OH. 45267-0536, USA.

IL-1 α has been found to be elevated in chronic wounds. Following treatment with a hydroactive dressing (Cutinova cavity), accelerated wound healing was associated with reduced IL-1 α levels. This mechanism is thought to be an IL-1 α mediated upregulation of matrix metalloproteinase-1, which, when present at elevated levels, prevents wound resolution¹³. There are conflicting reports of the role of VEGF in wound and bone healing literature. A systematic review of vascularity and healing of bone fractures concluded that VEGF was necessary for bone healing¹⁴. However, wound healing in diabetic mice can be improved by gene transfer of angiopoietin-1, without increased VEGF expression¹⁵ and blockade of VEGF has been proposed as a potential mechanism for improving cartilage healing in osteoarthritis¹⁶.

Analysis was done to compare the effects of InterX to the effects of TENS and demonstrated a statistically significant greater response in the InterX group, compared to the TENS group (Fig.3).



Conclusion:

This comparative study clearly demonstrates that both TENS and InterX influence similar responses in relation to lymphocyte metabolism and cytokine up and down regulation. However, the magnitude of that influence is significantly different with InterX having a much greater response in all parameters measured. TENS has a relatively small effect in some cases and a negligible effect in others.

The necessary level of lymphocyte activation needed to elicit a meaningful response to physiological stresses such as energy is unknown, though previous literature has demonstrated that TENS does not have an anti-inflammatory effect⁷. InterX has previously been shown to have an anti-inflammatory effect². The significant differences in the response to the two different types of stimulation may be explained by the fact that InterX delivers a much higher intensity and density of current than TENS as well as delivering stimulation specifically to optimal treatment points.

These data provide only a partial explanation of the inflammatory response and studies relating to circulation and lymphatic flow are warranted to further understand the complex interaction of mechanisms that can elicit a reduction of inflammation in patients. However, this study confirms the hypothesis that InterX elicits greater physiological response than TENS.

References:

1. Gorodetskyi I.G., et al. Non-invasive interactive neurostimulation in the post-operative recovery of patients with a trochanteric fracture of the femur. J Bone Joint Surg (Br) 2007;89-B:1488-94
2. Gorodetskyi I.G., et al: The effects of non-invasive interactive neurostimulation on pain and edema during post-surgical rehabilitation following internal fixation of unstable bi-malleolar ankle fractures. J Foot Ankle Surg 2010; 49:432-437
3. Maale G, Gamez M. The effects of a handheld, cutaneous, portable neurostimulator using two concentric electrodes with signals that are damped, bi-phasic oscillatory waveform which use the skin as a conduit in patients with chronic severe pain from large orthopaedic procedures. 18th Annual Symposium International Society for Technology in Anthroplasty, 2005:143
4. K. Gilbey; J. McKenna; E. Gardner; M.K. Hoeger Bement; M Johnson; K.A. Sluka; A.G. Leacock; D.M. Walsh. A Preliminary Investigation of the Hypoalgesic Effects of InterX using Experimental Pain: APTA CSM 2010, February 16-20, 2010, presented as a poster
5. Sluka K.A., Walsh D. Transcutaneous electrical nerve stimulation: basic science mechanisms and clinical effectiveness. J Pain 2003; 4:109-21
6. Han J S, Acupuncture: neuropeptide release produced by electrical stimulation of different frequencies. Trends in Neurosciences, Vol. 26, No.1, January 2003
7. Resende, M.A., et al, Local transcutaneous electrical stimulation (TENS) effects in experimental inflammatory edema and pain, European Journal of Pharmacology 504 (2004) 217– 222
8. Eizirik, D.L. et al, Human interleukin-1 beta stimulation of insulin release from rat pancreatic islets is accompanied by an increase of mitochondrial oxidative events. Diabetologia, 1989.32(11):p.769-73
9. New, K.J., et al, Effect of hypopolysaccharides and cytokines on oxidative metabolism in neonatal rat hepatocytes. J Pediatr. Surg., 2001 36(2): p338-40
10. Epstein, H., et al, innate immunity has a dual effect on vascular healing: supression and aggravation of neointimal formation and remodeling post-endotoxin challenge. Atherosclerosis, 2008. 199(1): p. 41-6
11. Yang, X., et al, Callus mineralization and maturation are delayed during fracture healing in interleukin-6 knockout mice. Bone 2007. 41(6): p. 928-36
12. Swartz, K.R., et al, interleukin-6 promotes post-traumatic healing in the central nervous system. Brain Res., 2001. 896(1-2): p.86-95
13. Barone, E.J.,et al, interleukin-1 alpha and collagenase activity are elevated in chronic wounds. Plast. Reconstr. surg. 1996. 102(4): p 1023-7
14. Keramaris, N.C., et al, Fracture vascularity and bone healing: a systematic review of the role of VEGF. Injury 2008. 39 (Suppl. 2): p. S45-57
15. Bitto, A., et al, Angiopoietin-1 gene transfer improves impaired wound healing in genetically diabetic mice without increasing VEGF expression. Clin. Sci. (lond), 2008. 114
16. Matsumoto, T., G.M. et al, blocking VEGF as a potential approach to improve cartilage healing after osteoarthritis. Musculoskelet. Neuronal Interact., 2008. 8(4): p:316-7