



INTERX5002

NRG
NEURO RESOURCE GROUP®

"InterX Technology is non-invasive, easy to use and gets results."

Dr. Gerhard Maale
Orthopedic Oncologist,
Presbyterian Hospital, Dallas

INTERX® 1000

Non-drug pain relief at your finger tips

Conveniently designed, safe and effective, the InterX 1000 can be used to supplement treatment with the InterX 5002 professional model, or as a home-use product. The InterX 5002 is a new option for the effective relief and management of a variety of painful injuries and chronic conditions.

INTERACTIVE TECHNOLOGY, PERSONAL APPLICATION

The InterX 1000 combines interactive technology with simple preset settings that allow immediate access to a range of treatment protocols.

- Five defined preset settings
- Fully portable
- Simple to use
- Convenient and safe
- Effective across a range of conditions



FLEXIBLE TREATMENT OPTIONS

Treatment with the InterX 1000 can be prescribed independently or in conjunction with other modalities. Typical treatment sessions range from 5 to 30 minutes.

PRODUCT SPECIFICATIONS

Size	220mm X 60mm X 45mm
Weight	Approx. 315 grams (11 ounces)
Power Source	4 AA DC Alkaline Batteries
Pulse Duration	10 – 500 µS
Pulse Frequency	15 – 350 Hz
Output Voltage	20 – 450 V
Output Current	2 – 90 mA
Waveform	pulsed, damped bi-phasic sinusoidal



PAIN MANAGEMENT SPECIALISTS HAVE TOLD US THAT THE INTERX 5002 ENHANCES REHABILITATION BY:

- Reducing acute and chronic pain
- Increasing post-injury range-of-motion
- Enhancing post surgical rehabilitation
- Improving muscle relaxation and reducing muscle cramps and spasms
- Enhancing neuromuscular re-education



Contact NRG directly or your local distributor to discover how your patients can benefit from non-drug, non-invasive pain relief with the InterX 5002.
Patent Pending

Email: info@nrg-unlimited.com US: 972 665 1810 UK: +44 (0)870 420 3204 www.nrg-unlimited.com



INTERX® 5002

*A dynamically different, clinical application
for the treatment of pain and injury*



"We perhaps have a singular opportunity to re-engineer post operative care – and I really mean that."

Dr. Arvind Peddada
Medical Director of Rehabilitation Services,
Presbyterian Hospital, Dallas

INTERACTIVE TECHNOLOGY

Designed specifically for pain and rehabilitation specialists, the InterX 5002 combines interactive technology with a new user interface that provides stimulation patterns to treat a full range of injuries and conditions.



The InterX 5002's unique interactive stimulation provides non-invasive, non-drug pain relief with simple clinical and post-surgical treatment applications.

INTELLIGENT SOFTWARE FOR OPTIMAL TREATMENT

The new graphical display provides users quick and easy access to treatment applications for acute and chronic conditions, as well as new treatment cycle programs.

The InterX 5002 has an additional "activity reading" feature that uses advanced numerical readings to objectively indicate optimal treatment points. The software measures and displays feedback from the impedance sensitive waveform and gives a numerical indication of where on the patient's body to treat—ensuring the most effective, targeted treatment.




OTHER

FEATURES INCLUDE:

- Fixed stainless steel electrodes that can be moved over the skin to identify optimal treatment points
- Impedance sensitive waveform that provides targeted stimulation
- Constantly changing impulses, avoiding accommodation
- New range of accessory electrodes
- 10-minute programs allow for optimal treatment time
- New unattended treatment options
- Athletic Injury Specific mode for sports specialists






"I have been using InterX for over a year and have been very impressed with the results. I hope the device becomes more widely available and other doctors and therapists benefit from using InterX as an aid in the management of chronic musculoskeletal pain."

Dr. Steve Motto
Sports and musculoskeletal physician,
London Bridge Hospital

SPECIFIC TREATMENT, EFFECTIVE RESULTS

THE INTERX 5002 AT A GLANCE

- Specifically designed for the management of chronic and acute pain, as well as rehabilitation from injury
- Different 'condition' stage applications provide easy access to preset stimulation patterns
- Fully portable, battery-operated device, weighing only 11 ounces (approx. 315 grams)
- Effective treatment times from 5 to 30 minutes
- Complete line of accessory electrodes to aid in therapy delivery
- Minimal set-up time, no skin preparation and no gel needed



The InterX 5002 is 510(k) cleared by the U.S. FDA for pain relief and pain management and carries the European CE mark.

"I was naturally skeptical, but I've used InterX for the past year and from the positive results I've seen, I now use it regularly as part of my treatment. InterX is being used more and more by my team and we have definitely been converted."

Kim Hobbs
Physiotherapist (PT),
Sports Dimensions Ltd., London

TARGETED APPLICATIONS FOR ADVANCED REHABILITATION

THE INTERX 5002 IS EFFECTIVE FOR:

- Reducing rehabilitation time from acute and chronic injuries
- The management of post surgical pain
- Pain relief and the resulting increase in range of motion
- Muscle relaxation, reducing muscle cramps and spasms
- Enhancing neuromuscular re-education

The InterX 5002 can be used effectively as an independent therapy or in conjunction with other modalities.

The InterX's ability to target specific areas of low impedance allows immediate access to highly effective treatment sites, often enhancing the rate of recovery.



EXPAND YOUR TREATMENT OPTIONS WITH THE NEW RANGE OF ACCESSORY ELECTRODES

- Tailored designs for a variety of applications, including unattended use
- Different sizes for specific surface areas
- Same interactive neurostimulation

CLINICALLY PROVEN*

- Significant decrease in average pain scores in acute and chronic patients
- Immediate increase in range-of-motion as a result of pain relief
- Reduction in estimated rehabilitation time frames



*Clinical trial data is available. Please contact NRG for details.